


The Establishment of Cannabinoid Family Receptors Screening Assays Facilitate Novel Drug Discovery

Xiaoyan Wang, dandan Hu, Tiejun (TJ) Bing, GPCR Department, ICE Bioscience Inc., Beijing, China
ICE Bioscience, INC. Building 14, Yard 18, Kechuang 13th Street, Beijing, 100176
Email: bingtj@ice-biosci.com




ICE Innovative GPCR Explorer
ICE Bioscience

Introduction

The CB1 receptor, a key component of the endocannabinoid system (ECS), is crucial in appetite control, pain perception, mood regulation, obesity, and diabetes. In oncology, the ECS plays a pivotal role in tumors' growth, development, and metastasis, with the activation or inhibition of the CB1 potentially exerting substantial impacts on cancer treatment and prevention. In vitro and in vivo modes are established to understand the signaling pathways of CB1, interactions with other cellular components and the development of more selective novel drugs with fewer side effects. CB1 receptor activation inhibits Forskolin-stimulated adenylate cyclase by G protein (Gai/o) and increases the phosphorylation of extracellular signal-regulated kinase 1/2 (pERK1/2) via G protein-dependent and β -arrestin-dependent pathways. Here we constructed an integrated experimental cascade, including in vitro HTRF cAMP assay and β -arrestin2 NanoBIT assay, to conduct the high throughput hit-to-lead compound screening, agonist, antagonist and inverse agonist validation. Meanwhile, several CDX modes are constructed for promisingly conducting in-vivo experiments and biomarker detection. Thus, ICE supports multiple approaches for helping the drug discovery and development of CB1 to facilitate the treatment of multiple diseases.

CB1 Screening Cascade



GPCR Overexpression Cell Line Construction




Figure 1. GPCR cell line construction by Flp-in system.




Figure 2. GPCR cell line construction by Lentivirus system.

Results

1. GPCR HTRF cAMP assay




Figure 3. HTRF cAMP assay was constructed for compound screening. The principle of this assay is grounded in HTRF (Homogeneous Time-Resolved Fluorescence) technology. It is formatted as an immunoassay for quantifying cAMP levels, which are generated in response to the modulation of adenylyl cyclase activity by GPCRs: Cannabinoid receptor CB1 and CB2 stable cell lines were constructed for agonist, inverse agonist and antagonist validation.

2. β -arrestin NanoBIT assay





Figure 4. Monitoring recruitment of the GPCR and β -arrestin in cells. A. The NanoBIT[®] assay operates on the principle that when two protein tags are in proximity to one another, the energy from a bioluminescent luciferase donor can be transferred non-radiatively to an acceptor tag, resulting in the emission of light from the acceptor. Transfer of energy is dependent on the distance and orientation of the tags to each other. Consequently, the assay can be effectively utilized to detect and measure the interactions occurring between proteins. B. The NanoBIT assay was used to characterize the CB1 and β -arrestin recruitment. C. C-C motif chemokine receptor 8 plays an important role in the migration and positioning of immune cells. The CCR8 and beta-arrestin1/2 recruitment was analyzed by NanoBIT assay. D. GPR35 is a multi-functional orphan receptor, a potential target for drug development and the treatment of inflammatory responses and cancer diseases. GPR35 and beta-arrestin1/2 NanoBIT assay was developed for drug discovery and compounds screening.

3. ICE_Obesity Panel Establishment





Figure 5. Obesity-Relevant Target Profiling: A. GPCR receptors related to obesity and metabolic disease were chosen to construct the ICE_Obesity Panel. B&C. Different receptors stable cell lines were established for application in multiple screening assays. The GIPR stable cell line was used for HTRF cAMP assay, GPR75 stable cell line was applied in reporter assay to validate its function. D. GPR75 SPR binding assay by Biacore. The signal changes during protein association and dissociation were recorded, and the affinity and kinetic parameters could be fitted.

4. DIO MODE in Mouse



Summary

We have successfully established a stable CB1 cell line that is versatile for various functional assays, ideal for the compound screening process aimed at identifying potential clinical candidates. Additionally, we have developed a β -arrestin2 NanoBIT assay tailored for the validation of oncology-related targets such as CB1, GPR35, and CCR8. This assay is designed to assess the activity of diverse ligands, encompassing agonists, antagonists, and biased agonists. Our CB1 screening cascade can provide comprehensive compound evaluation across in vitro and in vivo platforms, thus serving as an efficient screening platform for new drug discovery.

References

- [1] World Obesity Federation, World Obesity Atlas 2023. <https://data.worldobesity.org/publications/?cat=1>
- [2] Howlett AC, Abood ME. CB1 and CB2 Receptor Pharmacology. *Adv Pharmacol.* 2017;80:169-206. doi: 10.1016/bs.apha.2017.03.007. Epub 2017 Jun 12. PMID: 28826534; PMCID: PMC5812699.
- [3] Mangal N, Erridge S, Habib N, Sadanandam A, Reebye V, Sodergren MH. Cannabinoids in the landscape of cancer. *J Cancer Res Clin Oncol.* 2021 Sep;147(9):2507-2534. doi: 10.1007/s00432-021-03710-7. Epub 2021 Jul 14. PMID: 34259916; PMCID: PMC8310855.
- [4] Liu L, Doijen J, D'huys T, Verhaegen Y, Dehaen W, De Jonghe S, Schols D, Van Loy T. Biological characterization of ligands targeting the human CC chemokine receptor 8 (CCR8) reveals the biased signaling properties of small molecule agonists. *Biochem Pharmacol.* 2021 Jun;188:114565. doi: 10.1016/j.bcp.2021.114565. Epub 2021 Apr 17. PMID: 33872569.
- [5] Yansen Z, Lingang Z, Dali L, Mingyao L. Inflammatory bowel disease susceptible gene GPR35 promotes bowel inflammation in mice. *Yi Chuan.* 2021 Feb 16;43(2):169-181. doi: 10.16288/j.ycz.20-392. PMID: 33724219.
- [6] Zhang R, Xie X. Tools for GPCR drug discovery. *Acta Pharmacol Sin.* 2012 Mar;33(3):372-84. doi: 10.1038/aps.2011.173. Epub 2012 Jan 23. PMID: 22266728.