

A Biological and DMPK Integrated Platform for ADC Screening and Evaluation



Ying Meng, Li li, Xiao li, Min Kang, Lu Zhang, Jiabao Lv, Haiting Dai, Jiaying Zhang, and Tiejun Bing
International Department, ICE Bioscience InC.



Abstract

Antibody-drug conjugates (ADCs) are an innovative promising class of cancer therapeutics. ADCs are comprised three key components of an antibody (monoclonal antibody or bispecific antibody), a cytotoxic payload or other novel types of payload, and a linker. ADCs integrate chemotherapy and immunotherapy by combining the potency of payloads with the specificity of antibodies^[1]. ICE Bioscience has established a biological and DMPK integrated platform for ADC screening and evaluation. The platform is dedicated to support a comprehensive service portfolio with various aspects of ADC development projects.

Methods

1. Cytotoxic payload screening

- DNA Damage Repair inhibitor screening platform
- Topoisomerase I-mediated DNA relaxation assay

3. Antibody screening

- Target expression quantification
- Affinity: SPR
- Cell binding affinity
- Cell internalization

2. New MOA payload screening

- Targeted protein degradation assays
- STING agonist screening assays

4. ADC screening and evaluation

- Cytotoxicity
- Bystander effect
- Drug resistant cell line and CDX model
- DMPK

Case study: Cytotoxic payload screening DNA Damage Repair inhibitor screening platform




Figure 1. DNA Damage Repair (DDR) inhibitor screening platform in ICE. A. DDR inhibitor concept. B. Targets involved in ICE's DDR inhibitor screening platform. C. ATR inhibitor and PARP1 inhibitor screening assays.

Case study: Cytotoxic payload screening Topoisomerase I-mediated DNA relaxation assay




Figure 2. Inhibitory activity of DXd on Topoisomerase I (Top1). DXd prevented recombinant hTop1 from converting supercoiled DNA to relaxed closed circular DNA.

Case study: New MOA payload screening Targeted protein degradation assays




Figure 3. Targeted protein degradation assays for protein-degrader payload screening. A. Classical WB. B. Digital WB. C. IFA by High content imaging. D. In-cell WB assays.

Case study: New MOA payload screening STING agonist screening assays




Figure 4. STING agonist payload screening assays. A. SITNG HTRF binding assay. B. THP1 reporter assay. C. THP1-Dual reporter assay. D. Human PBMC cytokine release assay.

Case study: Antibody screening

Quantitative analysis of antigen membrane expression




Figure 5. Quantitative analysis of HER2 antigen membrane expression in different cancer cell lines with high/medium/low/negative HER2 expression. A. Histogram of Calibration Bead populations. B. Standard curve. C. Quantitative analysis of HER2 membrane expression.

Affinity - SPR




Figure 6. SPR analysis of Trastuzumab (Herceptin) binding affinity to HER2-ECD antigen and Cetuximab (Erbix) binding affinity to EGFR-ECD antigen using Biacore 8K.

Cell binding




Figure 7. Cell binding assays. A. Cell binding activity of Trastuzumab in different cancer cell lines with different HER2 expression levels. B. Cellular affinity comparison of DS-8201 and the corresponding antibody Trastuzumab on MCF-7 cell line.

Cell internalization




Figure 8. Antibody internalization evaluation. A. Antibody internalization assessed by flow cytometry. B. Kinetic monitoring of antibody internalization with Incucyte. C. Internalized antibody co-localization with lysosome, detected by High content imaging system.

Case study: ADC screening and evaluation

Cytotoxicity




Figure 9. Cell killing activity of DS-8201 and DXd in NCI-N87 cell line. A. Cell growing curves monitored by Incucyte. B. IC50 shifted ~10 times lower after conjugation to antibody.

In vitro bystander effect detection assays




Figure 10. In vitro bystander effect assay of DS-8201 and RC-48 in engineered breast cancer cell lines detected by flow cytometry. T-DM1 is negative control ADC for bystander effect assay. A. Flow cytometry results of SKBR3(GFP)&MDA-MB-468(RFP) co-culture after ADC treatment. B. GFP and RFP cell count and cell count ratio in the co-culture system.

Case study: Drug resistant cell line and CDX model




Figure 11. Development of drug resistance melanoma CDX model with A375 Dabrafenib-R cell line. A. Generation of A375-Dabrafenib R cell line. B and C. CDX tumor growth (B) and body weight (C) curves after inoculation.

Case study: DS-8201 Plasma Stability




Figure 12. In vitro stability of DS-8201 in plasma. A. DAR stability of DS-8201 in plasma. B. Concentration of free payload and total Antibody (TAB) in DS-8201 treated plasma. DS-8201 was incubated in plasma at 37 °C with different time points.

Conclusions

With 14 years of experience in providing early drug discovery services from target validation to pre-clinical candidate identification, ICE Bioscience has established a biological and DMPK integrated platform for ADC screening and evaluation. The platform is dedicated to support a comprehensive service portfolio with various aspects of ADC development projects, including antibody screening (target validation, affinity, binding, internalization), new payload (Topoisomerase I inhibitors, tubulin binders, etc.) or linker-payload screening (ADC characterization, bystander effect, etc.), efficacy and safety evaluation of the ADCs (gene editing and drug resistant cell line generation, cytotoxic effect, MOA, DMPK, etc.), in the process of ADC screening or PCC evaluation.

Acknowledgements

We thank ICE Bioscience for providing the research fund of this ADC biological and DMPK integrated platform.

References

[1] Tsuchikama K, Anami Y, Ha SYY, Yamazaki CM. Exploring the next generation of antibody-drug conjugates. Nat Rev Clin Oncol. 2024 Mar;21(3):203-223.