

VAV1-Induced Proximity: Insights from Molecule Glue Profiling

Xiaolan Su, Lili Chai, Tiejun Bing.  
ICE Bioscience, INC. Building 16, Yard 18, Kechuang 13th Street, Beijing, 100176  
Email: \*liy@ice-biosci.com



Background

Molecular glue degraders (MGDs) induce the proximity of target proteins to E3 ubiquitin ligases, leading to target protein ubiquitination and subsequent degradation by the proteasome. Unlike conventional inhibitors that target enzymatic activity, these molecules enable the degradation of proteins previously deemed "undruggable," thereby rendering them "druggable". With a first-in-class VAV1-directed MGDs (MRT6160) enters in clinical phase I, VAV1 degraders showed the remarkable potential in Immunology and inflammatory diseases such as rheumatoid arthritis and colitis. Current advancements highlight its ability to reduce proinflammatory cytokine production, inhibit pathogenic T cell polarization, and mitigate autoimmune responses.

Our integrated platform offers a comprehensive approach to exploring VAV1-targeting MGDs through a suite of state-of-the-art in vitro assays. These assays include the detection of binary and ternary complexes, intracellular protein interactions, and target degradation analysis using HiBiT assays. Additionally, the platform incorporates functional evaluations of T and B cell activities, cytokine profiling, and proteomics studies. Together, these components provide deep insights into the selectivity and potential off-target effects of these novel therapeutic agents.

Assay Summary Table

Binary/Ternary Complex Formation Assay		IC50 or Kd, nM		
HTRF biochemical assay	Binary	1217		
	Ternary	hVAV1	hVAV2	hVAV3
		18.35	>10,000	>10,000
		mouse/rat	dog	cyno
		18.04/14.28	14.08	17.89
Spectral shift biophysical assay		Binary: Kd=511 Ternary: Kd=16		
Cellular NanoBRET assay		IC50=18.45		
VAV1 degradation		DC50, nM		
Jurkat by HiBiT		8.98		
JESS analysis	Jurkat	T cell	B cell	PBMC
	3.79	0.65	0.94	0.99
Functional assays		IC50,nM		
T cells	Proliferation (Flow cytometry)	1.48		
	CD69 detection (Flow cytometry)	0.32		
	IL-2 secretion	0.35		
B cells	CD69 detection (Flow cytometry)	0.77		
	IL-6 secretion	0.75		
TH17 cell	IL-17A secretion by ELISA	0.08		

Summary

The common challenges of VAV1-like targets, such as complex ternary formation, stringent selectivity, and functional validation, are precisely the strengths of our platform. By addressing these commonalities through robust assay development and scalable workflows, we provide a universal solution for discovering and optimizing MGDs with broad therapeutic potential. Our platform doesn’ t just target VAV1—it enables researchers to uncover and exploit the shared principles of molecular glue activity, driving innovation across a wide range of targets.

References

[1] Application of induced proximity for therapeutic discovery. Cacace, Angela et al. Cell Chemical Biology, Volume 31, Issue 6, 1036 – 1038.

[2]. Neurath, M.F. and Berg, L.J. (2024) ‘Vav1 as a putative therapeutic target in autoimmune and chronic inflammatory diseases’, Trends in Immunology, 45(8), pp. 580–596. doi:10.1016/j.it.2024.06.004.

[3]. A vav1-directed molecular glue degrader, MRT-6160, reduces joint inflammation in a collagen-induced arthritis autoimmune disease model (no date) ACR Meeting Abstracts. Available at: <https://acrabstracts.org/abstract/a-vav1-directed-molecular-glue-degrader-mrt-6160-reduces-joint-inflammation-in-a-collagen-induced-arthritis-autoimmune-disease-model/> (Accessed: 22 January 2025).

Binary/Ternary Complex Formation Assay

To comprehensively evaluate binary and ternary complex formation in molecular glue studies, we have developed and applied multiple complementary assays:




Figure 2: Example data of CRBN&MGD binding A) and CRBN/VAV1 ternary complex formation B).

Degradation and Selectivity Analysis




Figure 5: A) and B) VAV1 degradation data in Jurkat-VAV1-HiBiT-KI cells. C) VAV1 degradation assay in Jurkat, T and B cells. D) Time-dependent VAV1 degradation: data from PBMC cells.

Immune Cell Activation Analysis

VAV1 associates with TCR and BCR signaling pathways and its degradation attenuates immune cell activation, proliferation, and cytokine production.




Figure 9: Example of VAV1 degrader that showed no off-target effects i iCESTP 90 Safety Panel