## Advancing HER2-Targeted ADC Drug Discovery through Resistant Cell Line **Generation and Cancer Cell Panel Analysis**

Lili Chai, Taotao Cao, Xue Yang, Wei Liu, Ying Meng and Tiejun Bing ICE Bioscience, INC. Building 14, Yard 18, Kechuang 13th Street, Beijing, China 100176 Email: bingtj@ice-biosci.com

#### Introduction

Antibody-drug conjugates (ADCs) represent a cutting-edge class of targeted cancer therapeutics, seamlessly integrating the high specificity of monoclonal antibodies with the cytotoxic power of chemotherapy agents. This strategic union facilitates the precise delivery of a potent cytotoxic payload to HER2-expressing cancer cells, where HER2 is a well-established therapeutic target. Despite the promising efficacy of HER2-directed ADCs in treating HER2-positive breast cancer, the emergence of resistance in tumor cells poses a significant obstacle to their therapeutic potential.



Figure 1. ICE database showed high expression level and Low expression level of HER2 in Cancer cell lines.



Figure 2. A, HER2 expression level in cancer cell lines detected by Flow. B, Binding assay in tumor cells with different HER2 expression levels. C, Comparison of Antibody and ADC binding EC50 in HER2+ SKOV3 cell. The result showed that NCI-N87, SKBR3, SKOV3 and HCC1954 cells have a high expression level of HER2 which makes them an ideal cell model for studies of breast cancer mono- and combination therapies.

## Payload Cytotoxicity and ADC Cell Panel Testing

Payload cytotoxicity test



Figure 3. Panels of 40-70 cancer cell lines were used to assess the cytotoxicity of Dxd and Exatecan. The results demonstrated that both compounds exhibited potent cytotoxicity across multiple cancer types.

#### **DS8201 efficacy in Cancer cell panel**



Figure 4. DS8201 showed antiproliferative activity across various cancer types, including breast, gastric, and non-small cell lung cancers. While effective in many HER2-overexpressing cell lines (e.g., AU565, HCC1954), some HER2-high-expressing lines exhibited weaker responses (e.g., SK-OV-3). This highlights DS8201's therapeutic potential and the variability in efficacy based on cellular context and HER2 expression levels.

SK-OV-3

#### **Tag-expression Reporter Tumor Cell Panel**

#### In vitro bystander effect assay

| Aa <sup>+</sup> cells          | В   |  |   |                  |  |  |
|--------------------------------|-----|--|---|------------------|--|--|
| Ag <sup>-</sup> cells          | 1   | Detection methods  | ADC bystander effect assays   | #                |  |  |
|                                |     | Flow cytometry   | Ag <sup>+</sup> /Ag <sup>-</sup> cell co-cultures with GFP/RFP  | 1                |  |  |
| SKBR3(                         |     | Flow cytometry   | Ag <sup>+</sup> /Ag <sup>-</sup> cell co-cultures without GFP/RFP   | 2                |  |  |
|                                |     | IncuCyte   | Ag <sup>+</sup> /Ag <sup>-</sup> cell co-cultures with GFP/RFP  | 3                |  |  |
| MDA-MB-4                       | đ   | IncuCyte (Realtime)  | Ag <sup>+</sup> /Ag <sup>-</sup> cell co-cultures with GFP/RFP  | 4                |  |  |
|                                | U   | IncuCyte and<br>Microplate Reader  | Co-cultures with Ag <sup>+</sup> -GFP/RFP cell<br>and Ag <sup>-</sup> Luc cell  | 5                |  |  |
| SKBR3(0<br>MDA-MB-4<br>co-cult |     | Microplate Reader  | Co-cultures with Ag <sup>+</sup> cell and Ag <sup>-</sup> Luc cell  | 6                |  |  |
|                                |     | Microplate Reader  | Medium transfer   | 7                |  |  |
| •                              | GFP | IncuCyte (Realtime)<br>IncuCyte and<br>Microplate Reader<br>Microplate Reader<br>Microplate Reader | Ag <sup>+</sup> /Ag <sup>-</sup> cell co-cultures with GFP/RFP<br>Co-cultures with Ag <sup>+</sup> -GFP/RFP cell<br>and Ag <sup>-</sup> Luc cell<br>Co-cultures with Ag <sup>+</sup> cell and Ag <sup>-</sup> Luc cell<br>Medium transfer | 4<br>5<br>6<br>7 |  |  |



Figure 5. The table sums up the methods to test the bystander effect, like flow cytometry, IncuCyte, and CTG. Figure B shows the bystander effect of DS8201 in co-cultures of HER2-positive SKBR3-GFP and HER2-negative MDA-MB-468-RFP cells detected by flow cytometry. Figure C gives the dose-curve effect of DS8201's bystander effect. Results show DS8201 has a strong bystander effect.

### Tag-expression cancer cell panel for bystander assays

| Cancer cell type  | GFP cells        | High expression targets                 | Cancer cell type        | Cell lines (ready)      |
|-------------------|------------------|---|-------------------------|-------------------------|
| Breast cancer     | SKBR3-GFP        | HER2, TROP2                             | Breast cancer           | MDA-MB-231-Luc          |
| Breast cancer     | HCC1954-GEP      | HER2 TROP2 MUC1 CDH3                    | Breast cancer           | MDA-MB-468-Luc          |
| Proact cancer     |                  |   | Breast cancer           | MCF-7-Luc-mEGFP         |
|                   |                  | MUC1, D7-114, CD115                     | Breast cancer           | H1299-Luc-mEGFP         |
| Breast cancer     | ZR-75-1-GEP      | MUC1, LIV-1, B/-H4                      | Brain Carcinoma         | LN299-Luc-mEGFP         |
| Ovarian cancer    | SKOV-3-GFP       | HER2, TROP2, Nectin4, B7-H4             | Brain Carcinoma         | U251-Luc-mEGFP          |
| Ovarian cancer    | OVCAR-3-GFP      | MSLN, B7-H4, MUC16, CLDN6               | Brain Carcinoma         | CT2A-Luc-mEGFP          |
| Gastric Cancer    | NCI-N87-GEP      | HER2 TROP2                              | Colorectal Carcinoma    | DLD-1-BRCA2-KO-Luc-mEGF |
|                   |                  |   | Colorectal Carcinoma    | HCT116-Luc-mEGFP        |
| Gastric Cancer    | NUGC-4-GFP       | ULDN 18.2, HER3                         | Lung cancer             | NCI-H1975-Luc-mEGFP     |
| Pancreatic cancer | BXPC-3-EGFP      | MET, CEACAM6, CEACAM5, CDCP1, EPHA2     | Lung cancer             | PC-9-Luc                |
|                   |                  |   | Lung cancer             | A549-Luc-mEGFP          |
|                   |                  |   | Lung cancer             | H1299-Luc-mEGFP         |
|                   |                  | Low expression terrets                  | Liver Carcinoma         | HepG2-Luc               |
| Cancer cell type  | RFP cells        | Low expression targets                  | Liver Carcinoma         | HepG2-Luc-mEGFP         |
| Breast cancer     | MDA-MB-468-RFP   | HER2, Nectin4, EGFR, B7-H4, NOTCH3      | Hepatoma                | HepG2-Luc               |
| Broast cancor     | MCE-7-REP        | FOLR1 Nectin4 FGER ROR1 PTK7            | Pancreatic cancer       | MIA-PaCa-2-Luc          |
|                   |                  |   | Prostate Adenocarcinoma | PC-3-Luc                |
| Lung cancer       | NCI-H1/92-RFP    | FOLR1, Nectin4, MUC1, ERBB3, CEACAM6    | Prostate Adenocarcinoma | PANC-1-Luc-mEGFP        |
| Lung cancer       | NCI_H226-RFP     | DLL3, TACSTD2, Nectin-4, CEACAM6, ITGB6 | Prostate Adenocarcinoma | BXPC-3-Luc-mEGFP        |
| Colon cancer      | SW620-RFP        | TACSTD2, CEACAM5                        | Prostate Adenocarcinoma | MIA PaCa-2-Luc-mEGFP    |
| Paparaatia aanaar | MIA-DoCo-2-DED   |   | Leukemia                | Jurkat-Luc              |
| Pancieatic cancer | IVIIA-FdGd-Z-KFF | GLUNTO.Z, TEGDU, 314                    | Leukemia                | K562-Luc                |
|                   |                  |   | Lymphoma                | Raji-Luc                |



| erative impact of DS8201 highlights its relative selectivity for HER2 overexpressing cell lines |   |                  |                      |                          |            |                    |  |  |  |  |  |
|---|---|------------------|----------------------|--------------------------|------------|--------------------|--|--|--|--|--|
| Cancer Type   | nTPM value for<br>ERBB2 (Data from<br>The Human Protein<br>Atlas) | HER2<br>Receptor | Estrogen<br>Receptor | Progesterone<br>Receptor | R_IC50, nM | Max<br>Inhibition% |  |  |  |  |  |
| Breast<br>adenocarcinoma  | 3987.3  | +++              | -                    | -                        | 0.21       | 95.98              |  |  |  |  |  |
| Breast<br>adenocarcinoma  | 801.3   | +                | -                    | -                        | 21.16      | 98.80              |  |  |  |  |  |
| Breast carcinoma  | 403.1   | +                | -                    | -                        | >100       | 3.20               |  |  |  |  |  |
| Breast carcinoma  | 3827.8  | +++              | -                    | -                        | 0.20       | 70.34              |  |  |  |  |  |
| Breast ductal carcinoma   | 3431.8  | +++              | -                    | -                        | 0.32       | 55.70              |  |  |  |  |  |
| Breast ductal carcinoma   | 1854.2  | ++               | -                    | -                        | 0.22       | 42.62              |  |  |  |  |  |
| Breast ductal carcinoma   | 3990.1  | +++              | -                    | -                        | 2.25       | 40.04              |  |  |  |  |  |
| Breast ductal carcinoma   | 2402.3  | +++              | -                    | -                        | 1.87       | 73.39              |  |  |  |  |  |
| Breast ductal<br>carcinoma  | 4058.3  | +++              | -                    | -                        | 0.15       | 75.49              |  |  |  |  |  |
| Breast metaplastic<br>carcinoma   | 2321.0  | +++              | -                    | -                        | 0.05       | 75.12              |  |  |  |  |  |
| Gastric tubular<br>adenocarcinoma   | 1288.3  | ++               | -                    | -                        | 0.45       | 53.86              |  |  |  |  |  |
| Non-small<br>cell lung cancer   | 1050.2  | ++               | -                    | -                        | 0.22       | 30.23              |  |  |  |  |  |
| Non-small<br>cell lung cancer   | 290.7   | +                | -                    | -                        | 1.55       | 85.08              |  |  |  |  |  |
| Non-small<br>cell lung cancer   | 3109.5  | +++              | -                    | -                        | 0.69       | 89.16              |  |  |  |  |  |
| Ovarian<br>adenocarcinoma   | 2047.9  | +++              | -                    | -                        | >100       | 15.50              |  |  |  |  |  |



Figure 6. Three ADC-related resistant cell lines for three cancer types (Colon, Gastric, and Ovarian) were generated via procedure A. DS8201 and three other payload molecules were evaluated in the resistant cell proliferation assays (B). Results indicated that DS8201 showed limited efficacy in the DLD-1 and DLD-1 Exatecan R cell lines, as well as in the SKOV3 and SKOV3 Dxd R cell lines. Conversely, MMAE exhibited minimal resistance across all three resistant cell lines compared to Dxd and Exatecan, suggesting distinct resistance mechanisms for these compounds.

In this study, we detected HER2 expression levels in various tumor cells using ICE's internal database and flow cytometry. We also screened for cancer types and cell lines sensitive to ADC payloads to establish ADC drug models. Additionally, tagged negative and positive cells were constructed to study ADC bystander effects, and resistant cell lines were developed from sensitive ones to investigate resistance and provide models to overcome it. Overall, the ICE platform offers comprehensive evaluation for ADC drug efficacy.

# Innovative CRO<sup>+</sup>Explorer Abstract Number: 1710

## ADC Resistant Cell Line

#### Summary