

STAT6: A Transcriptional Driver in Type 2 Immunity and a Rising Therapeutic Target

Tao Li, Qian Wang, Zhaoxia Yin, Yixiao Zhao, Qiang Xia, Xiaojian Wang, Tiejun Bing, Fujun Zhang, Ben Wei, Lili Chai | ICE Bioscience Inc., Beijing, China.

Email: bingtj@ice-biosci.com



Download PDF



Poster Number: P056

Abstract

Signal Transducer and Activator of Transcription 6 (STAT6) is a pivotal transcription factor activated downstream of the interleukin-4 (IL-4) and interleukin-13 (IL-13) signaling pathways. It governs the expression of key genes involved in type 2 immune responses, including Th2 cell differentiation, IgE class switching, and M2 macrophage polarization. Dysregulated STAT6 activity has been implicated in a range of diseases, from asthma and atopic dermatitis to fibrosis and immune-evasive tumors.

With the clinical advancement of STAT6-targeting drugs like Kymera's KT-621 and Recludix's REX-8756, the demand for precise STAT6 binding profiling is rising—supporting target validation, mechanism studies, and therapeutic development.

Representative STAT6 Inhibitors:


STAT6-IN-3: A small-molecule inhibitor blocking STAT6 phosphorylation or DNA binding for Th2-driven diseases.

AK-1690: A highly selective STAT6 PROTAC degrader with a DC50 of 1 nM, developed by the University of Michigan/ShaoMeng Wang group.

KT-621: An oral STAT6 PROTAC degrader with picomolar potency against IL-4/IL-13 signaling, developed by Kymera Therapeutics for atopic dermatitis and asthma.

I-2: The positive control in this article is molecule 2 from the Kymera Therapeutics patent (WO2025049820).



1 Binary/Ternary binder screening and orthogonal validation



Immobilized ligand	Analyte	1:1 binding ka (1/Ms)	kd (1/s)	KD (M)
STAT6	AK-1690	6.43E+05	2.53E-02	3.94E-08
STAT6	STAT6-IN-3	3.76E+06	1.30E-01	3.47E-08
STAT6	I-2	2.87E+06	2.56E-02	8.90E-09

Figure 1. Assay portfolio for STAT6 binder screening and ternary complex profiling. (A) HTRF assay for STAT6 binding screen. The IL-4Ra-probe can screen for small molecules and PROTACs that target the IL-4Ra binding site, such as STAT6-IN3 and AK1690. The M-probe can screen for small molecules and PROTACs that bind to the allosteric site, such as the KT-621 series. (B) HTRF ternary complex profiling for STAT6 PROTAC (hook effect found). (C) Spectrum Shift (SpS) assay for cost-effective HTS of STAT6 binary and ternary binding. (D) SPR as the gold standard for STAT6 binder screening.

2 Assessment of PROTAC-Mediated STAT6 Degradation




3 Functional Readouts in Primary Human Cells

Assay	IL-2	Dupilumab	YM-341619
IL-4 triggered pSTAT6 in CD3+ T cells of PBMC (FACS)	100.00	100.00	100.00
IL-4 triggered pSTAT6 in CD3+ T cells of human whole blood (FACS)	100.00	100.00	100.00
IL-4 induced Periostin release in Human Bronchial Smooth Muscle Cell	100.00	100.00	100.00
IL-13 induced Periostin release in Human Bronchial Smooth Muscle Cell	100.00	100.00	100.00
Reporter assay in HEK293-STAT6-Luc2p	100.00	100.00	100.00
IL-4 induced CCL17 release in human PBMC	100.00	100.00	100.00
IL-13 induced CCL17 release in human PBMC	100.00	100.00	100.00
IL-4 induced CD23 expression in B cells	100.00	100.00	100.00
IL-13 induced CD23 expression in B cells	100.00	100.00	100.00

Figure 3. Functional assays in primary human cells. (A) FACS-based analysis of STAT6 phosphorylation following stimulation in T cells. (B) Luciferase reporter assay measuring STAT6 transcriptional activity. (C) Periostin release from human bronchial smooth muscle cells (HBSMC) upon stimulation. (D) TARC secretion from PBMCs upon stimulation. (E) CD23 expression in B cells under different stimulation conditions.

4 Profiling STAT Family Selectivity and Cross-Species Specificity



Target	I-2	AK-1690
STAT1	/	11400
STAT2	/	1820
STAT3	/	2920
STAT4	/	2360
STAT5A	/	6020
STAT5B	/	186
STAT6	7.97	33.7

Figure 4. STAT family selectivity and cross-species specificity profiling. (A) HTRF evaluation of I-2-mediated STAT/CRBN ternary complex formation. (B) SPR analysis of STAT family selectivity. (C) Western blot-based assessment of STAT degradation selectivity. (D) Reporter assays across different STAT family members. (E) STAT phosphorylation measured by FACS, IFN- α -pSTAT1/2 and IL-12-pSTAT4 assays were ongoing. (F) Cross-species specificity evaluated by HTRF ternary complex formation, SPR binding, and degradation assays using several known STAT6 inhibitors/PROTACs.

5 Off-Target Profiling: Western, Proteomics, ICESTP Safety Panel & HiBit




Figure 5. Off-target profiling of STAT6 degraders (A) I-2 did not show off-target degradation of common CRBN substrates. (B) DIA-based proteomics profiling of I-2 off-targets. (C) ICE-STP safety panel revealed off-target effects of I-2 on LCK, A2a, PR, RARa and NR1/NR2B.

6 Developability: ADME, PK & HERG




Figure 6. Developability Profiling of STAT6 Degraders. (A) Pharmacokinetic profile of the STAT6 degrader following IV and PO administration. (B) HERG current measured by manual patch clamp. I-2 showed no significant inhibition of hERG current. (C) In vitro ADME properties of the STAT6 Degraders.

Summary

ICE Biosciences' discovery platform integrates HTS-compatible screening, orthogonal biochemical and cellular assays, and mechanistic studies. Supported by a broad suite of established assays and customizable modules, it generates high-quality activity and safety data to advance drug discovery pipelines. Incorporating human primary cell assays, PK/PD, and safety profiling ensures translational relevance, accelerating the development of next-generation STAT6 therapeutics with clinically meaningful insights.