

# Automated, Functional Dose-Response Screening for Comprehensive Pharmacological Safety Assessment

Zhiying Nie, Yanan Zhao, Meihui Shi, Ziyue Cui, Qiang Xia, Tiejun Bing, Fujun Zhang, Ben Wei, Lili Chai

Email: bingtj@ice-biosci.com

ICE Bioscience Inc., Beijing, China.



Download PDF



Poster Number: P033

## Introduction

In drug discovery, proactive toxicological risk assessment of candidate compounds is crucial for avoiding clinical-stage failure and post-market withdrawal. However, traditional preclinical safety evaluation strategies have significant limitations: in vivo animal studies exhibit low human relevance (translatability) due to significant species differences and face dual pressures of ethics and cost; while conventional in vitro binding assays are inadequate for predicting functional biological effects. Notably, approximately 75% of clinical adverse drug reactions (ADRs) originate from dose-dependent off-target effects, highlighting an urgent need for more accurate early-risk identification tools. To address these challenges, the IQ DruSafe consortium—a preeminent alliance of global pharmaceutical companies—has championed the enhancement of preclinical predictive power through the expansion and refinement of secondary pharmacology screening strategies. This study is aligned with this initiative and aims to investigate the implementation of an advanced in vitro secondary pharmacology screening system utilizing functional assay formats. By generating richer pharmacological information beyond mere binding affinity, this strategy seeks to enable a more accurate and earlier identification of potential safety liabilities, thereby providing highly translatable safety insights for the optimization of lead compounds.

The ICESTP SAFETYPANEL™ 77 Dose Response platform provides a comprehensive functional screening strategy that integrates single-point primary screening with quantitative dose-response (curve-based) profiling. This strategy is supported by a broad suite of mechanism-relevant functional assay technologies, including FLIPR calcium flux assay, HTRF, ADP-Glo, fluorescence polarization (FP), and other mechanism-relevant functional formats. The entire workflow exemplifies the convergence of functional biology and advanced automation. By accelerating the safety profiling process, minimizing manual error, and delivering mechanistically insightful, quantitative datasets (e.g., IC<sub>50</sub>/EC<sub>50</sub>), the platform enables more reliable optimization of candidate compounds early in the drug discovery pipeline.




Figure 1. Mechanisms of drug action

## Methods




Figure 2. This study adopted a high-throughput screening strategy integrating multiple functional assay technologies (including FLIPR, HTRF, QPatch, ADP-Glo, FP, etc.). The strategy included two modes: single-point primary screening and multi-concentration dose-response curve screening. To achieve high throughput and high precision, the platform integrated automated instruments (such as acoustic liquid handling systems and multi-mode plate readers) to perform curve-based screening tests. We applied this strategy to test the dose-response effects of various antidepressants, pergolide, and osimertinib on a panel covering 77 important safety-related targets.

## Overview of the ICESTP SAFETYPANEL™ 77






Figure 3. The ICESTP SAFETYPANEL™ 77 comprises 116 assays, integrating both single-point and dose-response screening methodologies within a functional assay framework.

## Results



## Functional Assays and Physiological ATP Enhance Prediction of Compound Safety Profiles




Figure 6. Pergolide with functional assay format ICESTP SAFETYPANEL™ 77 Dose Response. Using Pergolide as a case study, we systematically evaluated its off-target profile using the ICESTP SAFETYPANEL™ 77 Dose Response assay. As shown in Figure 6a, the heatmap visualization revealed distinct inhibition and activation patterns across multiple targets—particularly showing strong agonist activity at 5-HT<sub>1A</sub>, V<sub>1a</sub>, and α<sub>1A</sub> receptors. Subsequent dose-response characterization (Figure 6b) confirmed pergolide acts as a partial agonist at α<sub>1A</sub> (LogEC<sub>50</sub> = 1.945, Top = 46.19). Finally, we integrated these off-target effects into a mechanistic network (Figure 6c), linking its promiscuous receptor activity (e.g., 5-HT<sub>2B</sub>, α<sub>1A</sub>) to clinical adverse events such as valvulopathy and vasoconstriction. This approach provides actionable insights for predictive drug safety assessment.




Figure 7. Physiological ATP (1 mM) Enables Clinically Relevant Prediction of Osimertinib's Kinase Off-Target Profiles. A key advantage of the ICESTP SAFETYPANEL™ 77 Dose Response kinase assay is the uniform use of 1 mM ATP—a physiologically relevant concentration that mirrors the intracellular milieu—thereby significantly enhancing the clinical translatability of off-target predictions. While osimertinib potentially inhibits INSR and FAK under low, non-physiological ATP (Km), its efficacy is drastically reduced under 1 mM ATP (INSR IC<sub>50</sub> increases to 1397 nM; FAK inhibition is nearly abolished). This stark contrast demonstrates that only low IC<sub>50</sub> values obtained under physiological ATP conditions are clinically meaningful, preventing the overestimation of off-target risks that can occur in non-physiological assays.




Figure 9. Core Advantages of Functional Dose-Response Curve Analysis. This figure illustrates three pivotal advantages of employing functional dose-response curve analysis in ICESTP SAFETYPANEL™ 77 Dose Response: Accurately distinguishing full from partial agonism (Fig. A) by analyzing curve shape and maximal efficacy, providing critical insights into the mode of target engagement; Identifying non-linear responses (Fig. B), such as bell-shaped curves or aberrant plateaus, which may signal potential off-target effects or toxicity, thereby preventing false positives/negatives; Avoiding solubility-driven misreads (Fig. C), where precipitation of compounds at higher concentrations causes an anomalous decrease in our response. Curve analysis proficiently identifies these artifacts caused by physico-chemical properties rather than true pharmacological activity, ensuring accurate data interpretation.

## Conclusion and discussion

This study employed the high-throughput ICESTP SAFETYPANEL™ 77 to systematically profile the off-target activities of multiple drugs. The platform's automated functional analysis accurately delineated poly-pharmacological mechanisms underlying clinical safety profiles. Utilizing physiological conditions (e.g., 1 mM ATP for kinases) prevented risk overestimation and increased data reliability. Providing human-relevant safety insights with high predictive value. These results validate the platform's utility for mechanistically informed, early safety screening.

Looking ahead, high-throughput functional safety screening platforms represented by the ICESTP SAFETYPANEL™ 77 are poised to be deeply embedded into the entire drug discovery workflow, from target validation to clinical candidate selection. Their application can span multiple critical decision points: identifying and eliminating high-risk chemical scaffolds during the hit-to-lead stage; iteratively assessing the improvement in safety windows of different derivatives during lead optimization; and providing important mechanistic supplementary data for non-clinical safety evaluation during candidate drug nomination. By proactively and scalably applying such highly predictive in vitro safety profiling to R&D pipelines, drug discovery teams can identify potential safety risks at earlier stages and optimize resource allocation. This will significantly reduce late-stage development attrition due to safety issues and accelerate the development of new drugs with greater clinical advantages and improved safety profiles.

## Reference

- [1] Bowes, J. et al. Reducing safety-related drug attrition: the use of in vitro pharmacological profiling. *Nat. Rev. Drug Discov.* 11, 909–922 (2012)
- [2] Valentini, J. P. et al. In vitro secondary pharmacological profiling: an IQ-DruSafe industry survey on current practices. *J. Pharmacol. Toxicol. Methods* 93, 7–14 (2018).
- [3] Richard J. Brennan. et al. The state of the art in secondary pharmacology and its impact on the safety of new medicines. *Nat. Rev. Drug Discov.* 23, 525-545 (2024)
- [4] Stahl, S. M. *Stahl's Essential Psychopharmacology: Prescriber's Guide* (7th ed.). Cambridge University Press. (2021). (Page: Amitriptyline)
- [5] Anagha K, Shihabudheen P, Uvais NA. Side Effect Profiles of Selective Serotonin Reuptake Inhibitors: A Cross-Sectional Study in a Naturalistic Setting. *Prim Care Companion CNS Disord.* 2021;23(4)
- [6] Fava, M., Rush, A. J., Thase, M. E., Clayton, A., Stahl, S. M., Pradko, J. F., & Johnston, J. A. 15 years of clinical experience with bupropion HCl: from bupropion SR to bupropion XL. *Primary care companion to the Journal of clinical psychiatry*, 7(3), 106–113. (2005).
- [7] Maciag, M., & Karamyan, V. T. Enzymes in secondary pharmacology screening panels: is there room for improvement?. *Nature reviews. Drug discovery*, 24(6), 480–481. (2025).